

Uranotest[®]
Uriscreen[®]

For decision-making in suspected cases
of urinary tract infections




urano[®]
vet

Uranotest® Uriscreen®

The quick, reliable, low-cost tool for correct decision-making in urinary tract infections.

Uranotest® Uriscreen® is a test for detecting bacteria and somatic cells in urine that makes it possible to diagnose urinary tract infections caused by catalase-positive bacteria (the main cause of urinary tract infections) and the presence of somatic cells in just 2 minutes.



In cases of suspected urinary tract infections, Uranotest® Uriscreen® enables us to decide which samples need to be sent to the laboratory for culture.


The main advantages of using Uranotest® Uriscreen®

- High sensitivity even with low bacterial counts as it detects concentrations from just 50,000 CFU/ml.
- Negative (-) predictive value greater than 95% (percentage of negatives with Uranotest® Uriscreen® that truly are negative).
- Result in less than 2 minutes.
- Economic, which means it can be used in screening for decision-making.
- No equipment needed.

Procedure

1

With the help of the pipette supplied, add a volume of urine to the tube so the level is between the minimum and maximum lines marked.



2

Add 4 drops of reagent.




3

Mix gently without stirring for 5 seconds until the urine turns blue.



4

Observe the formation of foam on the surface of the liquid.



The formation of some bubbles due to the shaking must not be considered a positive result. In any case, a consistent foam must form.

The detection of catalase-positive bacteria as an indicator of urinary tract infection

Catalase action (bacteria)	Expected result with Ura-notest® Uriscreeen®	Visual result	Interpretation	Recommendation
Catalase + <i>Staphylococcus pseudintermedius</i> <i>Corynebacterium diphtheriae</i> <i>Enterobacteriaceae</i> <i>Escherichia coli</i> <i>Citrobacter spp.</i> <i>Enterobacter spp.</i> <i>Klebsiella spp.</i> <i>Proteus spp.</i> <i>Salmonella spp.</i> <i>Pseudomonas spp.</i>	Positive	Foam is generated on the surface until a complete ring is formed.	There may be bacteria and/or leukocytes in the urine suggestive of urinary tract infection due to catalase-positive bacteria.	Send the sample to the laboratory for culture and antibiogram. Start treatment with a broad spectrum antibiotic effective for urinary tract infections (amoxi + clavulanic acid or cephalosporin). When you have the culture results, change the antibiotic for a more sensitive one or maintain if the current one is sufficiently sensitive.
Catalase – <i>Streptococcus spp.</i> <i>Enterococcus spp.</i>	Negative, although if the sample has more than 10 leucocytes per field as a consequence of the infection, the result may also be positive.	Foam does not form on the surface or an incomplete ring of foam forms.	Catalase-positive bacteria not present in the urine.	There is no need to send the sample to the Laboratory for a culture, as the test has a negative predictive value of 95 %.

Comparison of sensitivity and specificity of Uranotest® Uriscreen® and study of sediment under microscope verses urine culture

Screening n= 165 dogs and cats	Sensitivity	Specificity
Uranotest® Uriscreen®	89 %	71 %
Determination under microscope by expert staff of the presence of abnormal urinary sediment (pyuria, bacteriuria or both))	78 %	90 %

Detection limit	50.000 CFU/ml
	10 leucocytes/field

Causes of false positives and false negatives

False negatives	<p>Catalase-negative bacteria (Streptococcus, Enterococcus)</p> <p>Bacterial concentrations below the detection limit.</p>
False positives	<p>Bacteria that not does not grow in cultures (Chlamydia, Trichomonas)</p> <p>Haematuria</p> <p>Presence of cells shed from the bladder or kidneys in urine</p> <p>Presence of tumour cells</p>

BIBLIOGRAPHY

- Vail DM, Allen TA, Weiser G. Applicability of leukocyte esterase test strip in detection of canine pyuria. Journal of the American Veterinary Medical Association 1986; 189(11):1451-3.
- Weese SJ, Blondeau JM, Boothe D, et al. Antimicrobial use guidelines for treatment of urinary tract disease in dogs and cats: antimicrobial guidelines working group of the international society for companion animal infectious diseases. Veterinary Medicine International. 2011 article ID 263768.
- Bartges JW. Diagnosis of urinary tract infections. Veterinary Clinics of North America Small Animal Practice 2004; 34:923-933.
- Forester SD, Troy GC, Dalton MN, et al. Retrospective evaluation of urinary tract infection in 42 dogs with hyperadrenocorticism or diabetes mellitus or both. Journal of Veterinary Internal Medicine 1999; 13:557-560.
- Litster A, Moss S, Platell J, et al. Occult bacterial lower urinary tract infections in cats: Urinalysis and culture findings. Veterinary Microbiology 2009; 136:130-134.
- Lekcharoensuk C, Osborne CA, Lulich JP. Epidemiologic study of risk factors for lower urinary tract diseases in cats. Journal of American Veterinary Medical Association 2001; 218(9):1429-35.
- Peterson AL, Torres SM, Rendahl A, Kock SN. Frequency of urinary tract infection in dogs with inflammatory skin disorders treated with ciclosporin alone or in combination with glucocorticoid therapy: a retrospective study. Veterinary Dermatology 2012 Jun; 23(3):201-e43.
- Ling GV, Biberstein EL, Hirsh DC. Bacterial pathogens associated with urinary tract infections. Veterinary Clinics of North America Small Animal Practice 1980; 9:617-630.
- Schonbaum GR, Chance B. in The Enzymes, 2nd edition, P.D. Boyer, ed., Vol 13, pp. 363—408 (1976), Academic Press, New York.
- Nadler V, Goldberg I, Hochman A. Comparative study of bacterial catalases, Biochim. et Biophys. Acta 882: 234-241 (1986).
- Jarvinen KAJ, Determination of peroxidase in urine. British Medical Journal, 1: 379 (1958).
- Braude AI, Pittsburgh HB. Detection of urinary catalase by disk flotation, Journal Laboratory and Clinical Medicine. 57: 490-494 (1961).
- Lieberman I, Ove T. Enzyme activity levels in mammalian cell cultures, Journal of Biological Chemistry 223: 634-636 (1958).
- Braude AI. Catalase activity of infected urine, J. Clinical Investigation 38: 990 (1959).
- Millar L, Debuque L, Leialoha C, et al. Rapid enzymatic urine screening test to detect bacteriuria in pregnancy. Obstetrics and Gynecology 2000;95:601-604.
- Kvitko-White HL, Cook AK, Nabity MB, Zhang S, Lawhon SD. Evaluation of a catalase-based urine test for the detection of urinary tract infection in dogs and cats. Journal of Veterinary Internal Medicine 2013; 27:1379-1384.



Uranovet SL
Avda. Santa Eulalia, 2
Les Franqueses I Barcelona
www.uranovet.com