

Uranotest Dermatophytes



Dermatophytes

- Affect keratinised tissues: skin, nails and corneal stratum.
- Their diagnosis is very important due to their zoonotic potential.
- Prevalence is greatly increased due to the importing of puppies in a dubious state of health, an increase in adoptions and a rise in the habit of keeping rodents as pets

Presentation:

Box of 4 tests

Uranotest Dermatophytes: a new faster, more convenient and more reliable culture medium for the diagnosis of Dermatophytosis

Fewer false positives
due to growth of
saprophyte flora

Much simpler
seeding than with
the tube system

5 cm diameter dish

Easy collection of
the colonies for
re-seeding and/or
observation under
the microscope

Colonies grow
further apart and are
easier to identify

Quick colour
change (after as
little as 2 days)



Includes contrast medium for easy identification and differentiation of the species when observing them under the microscope

Species most frequently involved

Cats	Microsporum canis
Dogs	Microsporum canis Microsporum gypseum Trichophyton mentagrophytes
Rodents	Microsporum canis Trichophyton mentagrophytes

Tips for performing the test properly

Taking samples

- Washing the skin is only recommended in the case of contamination and the presence of scabs. If done, always use a non-medicated soap that does not have antifungal properties and dry well before taking the sample.
- With the help of a scalpel or forceps, select hair and desquamation both from the periphery and the centre of the lesion. Brittle and broken hairs, as well as those which display fluorescence when observed with Wood's Lamp, are the best samples.
- Take only a small amount of hair and scales, since excess can induce the growth of saprophytic flora.
- In asymptomatic cats, it is advisable to obtain a sample of epithelial desquamation of the suspicious area using a sterile toothbrush. Once the area has been brushed, the tips of the brush fibres can be cut and incubated directly onto the dish.

Sowing and incubation

- Sow onto the culture medium, ensuring good contact with its surface, but without "burying" the sample in the medium.
- The ideal incubation temperature is 27 - 28° C. If a incubation oven is not used, keep in a warm place away from the light. Very low or very high incubation temperatures may lead to false positive results.
- The lid is designed with 3 fins that allow air to enter the medium.
- The dish should be incubated preferably in a dark place with the lid facing downwards.

Interpreting the results

Change of colour	Time period	Colour of the colonies	Interpreting the results
None	After 12 days	No colonies	NEGATIVE
Yellow to red	Between 2 and 12 days	White colonies	*POSITIVE
Yellow to red	After 12 days	Brown, grey or greenish colonies	NEGATIVE. The change of colour is due to the growth of saprophyte flora which occurs after the recommended maximum reading time of 12 days
None	After 2 days	Brown, grey or greenish colonies	NEGATIVE

*Interpreting the results: after the 2nd day following seeding



Microsporum
Flat, white velvety colony.



Trichophyton
A dusty look colony, flatter than that which M. canis presents.

The culture medium must be observed daily for a maximum of 12 days to check the colour change and the growth of colonies.

Identification under the microscope

Although in most cases a positive result already makes it possible to prescribe an appropriate treatment, the Uranotest Dermatophytes dish has been designed to allow the colonies to be collected using an adhesive tape, so that they can subsequently be observed under a microscope and thus identify the species involved.

When the colony is 10 - 12 days old, the clear adhesive tape is laid over the colonies, pressing gently. The sample is transferred to a slide, on which a drop of contrast solution (supplied with the kit) has been previously placed. The adhesive tape itself acts as a cover slip.

Under the microscope, the macroconidias can be observed and identified by performing a differential diagnosis between the different species.

Microsporum canis

Abundant macroconidias in the shape of the keel of boat usually present over six septa. The microconidia are scarce, small and club-shaped.




Microsporum gypseum

Abundant thin-walled, elliptical macroconidias showing less than 6 septa. If they are present, the number of microconidia is very low, club-shaped and with a smooth wall.



Trichophyton mentagrophytes

Cigar-shaped macroconidias. Rounded microconidia grouped into a cluster.



Urano Vet SL

Parque Empresarial Mas Blau I
Calle Garrotxa, 10 -12 . Planta 1B
08820 El Prat de Llobregat I Barcelona
Phone: + 34 93 193 52 55
E info@uranovet.com
W www.uranovet.com