Uranotest Copro is a device that simplifies coprological examination, from the process of sample collection by the owner until the completion of the parasitological examination at the clinic.

Uranotest[®]Copro All in the same device

- 1 Collection
- 2 Dilution
- 3 Preservation
- 4 Filtration
- **5** Concentration

Without odours or handling of stools



Bottle containing Greenfix[®], a preservative that is non-toxic and non-aggressive to the environment, and that ensures the sample is preserved and kept in good condition for up to 15 days.

The same bottle incorporates a 260 micron filter that lets through eggs, larvae and cysts, while retaining large particles, thereby obtaining a clean sample for microscopic analysis.